

Photovolt Aquatest 1010

Karl Fischer Coulometric Moisture Titration System

- Extremely small footprint for laboratory bench, 5 inch width and lightweight, ~5 lbs.
- CE compliant, allowing easy entry into international markets
- Cabled keypad/display can be placed on bench, on unit or mounted on wall.
- RS-232 interface for optional printer or computer connection
- Totally greaseless fittings for glassware.
- Quick accurate moisture analysis, provides maximum isolation of solutions from atmospheric moisture.
- Photovolt's Karl-Fischer reagents available in kit and bulk packaging.

Applications

Pharmaceutical
Biomedical
Petrochemical

Chemical
Plastics
Utilities

Food & Beverage
Transportation
Aerospace

