

Photovolt Aquatest 2010

Karl Fischer Coulometric Moisture Titration System

- Full-featured, user friendly operating system, multiple RS-232C ports for system expansion, simultaneously interfacing both a balance and a computer.
- Quick, accurate moisture analysis, provides maximum isolation of solutions from atmospheric moisture.
- Optional Vaporizer analyzes challenging materials like lyophilized proteins, pharmaceutical powders, or biomedical, pharmaceutical and plastic industry solids.
- Photovolt's Karl-Fischer reagents, available in kit and bulk packaging.
- GLP software, includes statistics, diagnostics and reagent monitoring packages.
- CE-Mark compliant, allowing easy entry into international markets.

Applications

Pharmaceutical	Chemical	Food & Beverage
Biomedical	Plastics	Transportation
Petrochemical	Utilities	Aerospace

Photovolt Instruments

